

SPECTRUM DPL PLUS



2024 NEWEST DELICATE PULSED LIGHT

DPL Hair Removal + Filter DPL Skin Care



Create a new

skin rejuvenation revolution

DPL HR HANDLE

Spot size: 10*50mm

Wavelength: 650-950nm

Lamp lifespan: 1 million shots



FILTER DPL HANDLE

Spot size: 10*50mm






Wavelength: 420nm;530nm; 515 nm, 560 nm; 590 nm; 615 nm;
640 nm; 695 nm; 755 nm;

Lamp lifetime: 1 milllion shots




Handle 1: DPL HR

DPL Hair Removal is suitable for both men and women and can be performed on all skin types.




~~Handle 1: DPL HR~~

Difference between IPL & DPL hair removal



- Removes all the gentle and small hairs
- Uses only pulses of light that to in a controlled manner which essentially treats the small hairs too.
- Treats all types of hairs



- Remove only coarse hairs first and then a small number of tiny hairs.
- Continuous exposure to the light which is time-consuming and painful too.
- Can be used only on specific skin tones and hair colors

~~Handle 1: DPL HR~~



• 8HZ Fast Hair Removal

Handle 2: Filter DPL

With the Filter DPL, The Skin Rejuvenated is Easy.



Acne



Epidermis pigmentation



Epidermis pigmentation,
Rosacea and flushing



Pigmentation,
Vascularisation and skin care.



Dermis larger vessels,
Skin care,
Pigmentation and vascularization



Larger vessels in the face,
General ipl skin care,
Pigmentation and vascularization



Oil control,
Fine lines,
Skin texture improvement




Hair removal



Deep anti-aging

Principle

The longer the wavelength, the deeper the penetration



- 01** Fine lines and wrinkles 
- 02** Pigment spots, sun damage 
- 03** Enlarged pores 
- 04** Redness, couperose, rosacea 
- 05** Skin rejuvenation and acne 

~~Cooling Handle~~

- 4 - 20°C
Cooling



Spot Size Adaptors


Quick-Change Spot Size Adaptors

Quickly and easily change spot sizes with magnetic gold pieces that instantly adjust for comfortable treatment of both small and highly contoured areas.



OPT and MSP Technology


OPT - Optimal Pulse Technology



The OPT advantages:

- High peak power, shorter pulses - optimal for DPL skin treatments using photorejuvenation and the treatment of benign pigmented lesions.
- Homogeneous power over long pulse durations.

MSP - Multiple-Sequential Pulsing



MSP Technology:

- Multiple-Sequential Pulsing allows for cooling the skin in between pulses, reducing the chance of skin damage during treatments.


Advantages



Treats Acne, Rosacea, Pigmented and Vascular Lesions.
Also reduce age spots through photorejuvenation.




9 Filters, changeable in seconds.



Optimal Pulse Technology (OPT)
controls pulse shape and length.



-4-20°C for patient comfort
through continuous contact cooling.



Multiple Sequential Pulsing (MSP)
enables a safer procedure while cooling
between pulses.



DPL fast hair removal

UI Intertace



SDPL for fast DPL hair removal



FILTER DPL with 9 filters



DPL for professional hair removal



Parameter

Specification of DPL Plus

WAVE LENGTH
420-
1200nm(Filters)
560-950nm(HR)

FILTERS
420nm;530nm; 515 nm, 560 nm; 590 nm; 615
nm;
640 nm; 695 nm; 755 nm;

FLUENCE
5-25J/cm²

INPUT POWER
4000W

POWER SOURCE
Pure Sapphire

LAMP
Germany Imported lamp

SPOT SIZE
15*50mm

CONTACT COOLING TEMPERATURE
-4°C-20°Cadjustable

PACKAGE SIZE
65*55*121cm


COOLING SYSTEM
Built in water cooling +semi conductor
cooling +air cooling

WEIGHT
87kg

LCD SCREEN
10.4 inch true color touch screen control



~~Hair removal~~



Before



After Immediate effect

•We're confident in our results, so you can be too

Hair removal



Before



After Immediate effect

•We're confident in our results, so you can be too

Shrink pores



Before

After Immediate effect

•We're confident in our results, so you can be too

Rosacea



Before



After 2 treatments

•We're confident in our results, so you can be too

~~Freckle~~



Before



After 2 weeks

•We're confident in our results, so you can be too

~~Vascular~~



Before



After 4 weeks

•We're confident in our results, so you can be too

Acne inflammation



Before



After 2 treatments

•We're confident in our results, so you can be too

lass lip



Before

After Immediate effect

•We're confident in our results, so you can be too



REACH OUT TO OUR DEDICATED TEAM FOR ANY INQUIRIES, ASSISTANCE, OR INFORMATION YOU NEED.

 ADMIN@SKY-DERM.COM

 **UAE BRANCH:**
DUBAI: OFFICE NO. G01-606, DEEB INVESTMENT COMPLEX

 +971 56 256 5660 +971 504807275

 **EGYPT BRANCH:**
CAIRO: 40 ABDALLAH EL ARABY, TAYARAN
ST, NASR CITY, EGYPT

 +20 102 623 9460

 WWW.SKY-DERM.COM